

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E0775724B

Time Processed :8:56:09

Month-Day-Year Processed :03/02/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.573	8802	532	0.071
2	12.071	59698	2077	0.482
3	17.611	23002	3301	0.186
4	17.939	18306	1337	0.148
5	18.738	196365	21507	1.585
6	19.470	156831	15969	1.266
7	20.168	11914446	1141705	96.196
8	21.953	8206	609	0.066
Total		12385655	1187036	100.000