

Sample Name :Melittin,honey bee
Cat No. :RP20415
Lot No. : E0838325K
Time Processed :21:38:23
Month-Day-Year Processed :11/22/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.433	61420	7363	0.760
2	20.392	37799	3984	0.468
3	20.575	77688	8198	0.961
4	20.642	45473	7518	0.562
5	20.992	7800225	871142	96.479
6	21.375	48352	6361	0.598
7	21.567	13945	2345	0.172
Total		8084901	906911	100.000