


Mass Spectrum



Sample Information

Month-Day Processed : 10/14/24
Time Processed : 18:35:56
Injection Volume : 0.3
Sample Name : TAT peptide

Cat No. : RP20256
Lot No. : E0838424J
Theoretical MW: 1621.91

Observed MW : 1621.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH