


Sample Name :KRAS G12V Peptide (VVVGAVGVGK)
Cat No. : RP30229
Lot No. : E0871425F
Time Processed :10:24:37
Month-Day-Year Processed :06/08/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ20010139

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.658	25615	3228	0.348
2	9.767	20467	4936	0.278
3	9.933	6998684	1118804	95.062
4	10.192	45620	6615	0.620
5	10.300	27026	4896	0.367
6	10.675	29649	4464	0.403
7	11.183	60281	11300	0.819
8	11.433	19028	3484	0.258
9	15.233	114248	10766	1.552
10	24.858	21617	2605	0.294
Total		7362233	1171098	100.000