


Mass Spectrum



Sample Information		Equipment : ZI21010035	
Month-Day Processed :	07/10/23	Interface Bias :	+4.5 kV
Time Processed :	15:07:49	Drying Gas Flow :	5 L/min
Injection Volume :	0.4	T.Flow :	0.2 ml/min
Sample Name :	DYKDDDDK Peptide	B.conc :	50%H2O/50%MeOH
Cat No. :	RP10586		
Lot No. :	E0912623G		
Theoretical MW: 1012.97			
Observed MW : 1013.0			