

Sample Name:Hepatitis B Virus Core (128-140)

Cat No. : RP17313

Lot No. : E0961025F

Time Processed :14:27:49

Month-Day-Year Processed :07/09/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.138	24203	2415	0.090
2	12.287	10107	1837	0.037
3	12.410	32123	4874	0.119
4	13.355	26379382	1515073	97.801
5	13.772	444560	54823	1.648
6	14.002	82114	11003	0.304
Total		26972489	1590026	100.000