


Mass Spectrum



Sample Information

Month-Day Processed :	07/09/25	Equipment	: ZJ21010035
Time Processed :	13:44:56	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	Hepatitis B Virus Core (128-140)	T.Flow	: 0.2 ml/min
Cat No. :	RP17313	B.conc	: 50%H2O/50%MeOH
Lot No. :	E0961025F		
Theoretical MW:	1406.63		
Observed MW :	1406.6		