

Sample Name:OVA Peptide (323-339)
Cat No. : RP10610
Lot No. : E1059126B
Time Processed :23:56:02
Month-Day-Year Processed :02/12/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ25010054

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.050	10517615	1474262	98.866
2	8.433	39355	4711	0.370
3	8.558	20922	3695	0.197
4	8.792	60361	9744	0.567
Total		10638253	1492413	100.000