


Mass Spectrum



Sample Information

Month-Day Processed : 03/01/26
Time Processed : 8:05:00
Injection Volume : 0.3
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E1060826B
Theoretical MW : 4186.66
Observed MW : 4186.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment
Interface Bias
Drying Gas Flow
T.Flow
B.conc

: ZJ22010150
: +4.5 kV
: 5 L/min
: 0.2 ml/min
: 50%H2O/50%MeOH