

Sample Name :TAT peptide
Cat No. : RP20256
Lot No. : E1111725B
Time Processed :1:27:11
Month-Day-Year Processed :03/03/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ20010139

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	4.158	65633	10522	0.739
2	5.600	14762	2050	0.166
3	5.958	8644692	899562	97.278
4	6.158	30910	4702	0.348
5	7.033	65992	2431	0.743
6	8.067	39805	1683	0.448
7	8.283	24768	2856	0.279
Total		8886562	923807	100.000