


Mass Spectrum



Sample Information

Month-Day Processed : 03/03/25
Time Processed : 12:20:08
Injection Volume : 0.3
Sample Name : TAT peptide

Cat No. : RP20256
Lot No. : E111725B
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

:ESI

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min

Equipment

: ZJ21010035
B.conc : 50%H2O/50%MeOH