

Sample Name : Melittin,honey bee
Cat No. : RP20415
Lot No. : E1130922H
Time Processed :15:35:45
Month-Day-Year Processed :08/18/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
35.00	Pumps	B.Conc	95
37.01	Pumps	B.Conc	5
37.02	Pumps	B.Conc	5
37.20	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.994	24096	3852	0.252
2	20.521	60529	6116	0.632
3	20.636	121302	9321	1.267
4	21.072	9147281	1040565	95.548
5	21.616	134984	9960	1.410
6	21.820	21678	3134	0.226
7	22.031	63596	9165	0.664
Total		9573465	1082113	100.000