


Mass Spectrum



Sample Information

Month-Day Processed :	08/15/22	Equipment	:ZJ22010150
Time Processed :	2:28:23 PM	Interface Bias	:+4.5 kV
Injection Volume :	0.1	Drying Gas Flow	:5 L/min
Sample Name :	Melittin,honey bee	T.Flow	:0.2 ml/min
Cat No. :	RP20415	B.conc	:50%H2O/50%MeOH
Lot No. :	E1130922H		
Theoretical MW :	2847.50		
Observed MW :	2848.0		