


Mass Spectrum



Sample Information		Interface	:ESI	Equipment	:ZJ22010150
Month-Day Processed :	04/11/24	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Time Processed :	10:58:04	CDL Temp	:250	Drying Gas Flow	:5 L/min
Injection Volume :	0.2	Block Temp	:200	T.Flow	:0.2 ml/min
Sample Name :		[Gly22]-β-Amyloid (1-42), Arctic Mutation			
Cat No. :	RP30247CN	B.conc			
Lot No. :	E1139424D				
Theoretical MW :	4442.00				
Observed MW :	4442.0				