

Sample Name :[Gly22]- β -Amyloid (1-42), Arctic Mutation

Cat No. : RP30247CN

Lot No. : E1140024D

Time Processed :18:10:50

Month-Day-Year Processed :04/15/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.958	47480	1439	0.670
2	15.667	37563	1204	0.530
3	17.225	6908846	204961	97.564
4	18.408	87469	2529	1.235
Total		7081359	210132	100.000