

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E1227923J
Time Processed :20:24:21
Month-Day-Year Processed :11/02/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010014

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	13.817	17293	2322	0.677
2	14.392	2477335	280758	96.985
3	14.808	26754	1680	1.047
4	15.225	32960	3646	1.290
Total		2554341	288407	100.000