

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E1234323J

Time Processed :15:58:20

Month-Day-Year Processed :11/08/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ20010139

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	14.883	34731	4056	0.471
2	15.000	32264	4626	0.437
3	15.225	81568	9704	1.106
4	15.542	7071920	918862	95.876
5	15.967	123141	14577	1.669
6	16.433	5140	341	0.070
7	17.233	22184	2419	0.301
8	17.475	5145	797	0.070
Total		7376093	955381	100.000