


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	11/02/23	Interface	: ZJ21010035
Time Processed :	12:27:25	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	T.Flow	: 0.2 ml/min
Cat No. :	RP10769	B.conc	: 50% H ₂ O/50% MeOH
Lot No. :	E1234323J		
Theoretical MW :	3766.20		
Observed MW :	3765.9		