


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	11/08/23	Nebulizing Gas Flow :	1.5L/min	Interface Bias	:ZJ22010150
Time Processed :	10:05:26 AM	CDL Temp :	250	Drying Gas Flow	: +4.5 kV
Injection Volume :	0.1	Block Temp :	200	T.Flow	:5 L/min
Sample Name :	Amylin, human, amide			B.conc	:0.2 ml/min
Cat No. :	RP11278				:50% H2O/50% MeOH
Lot No. :	E1236623J				
Theoretical MW :	3903.28				
Observed MW :	3903.2				