


Mass Spectrum



Sample Information

Month-Day Processed :	02/27/24	Equipment	:ZJ22010150
Time Processed :	20:22:37	Interface Bias	:+4.5 kV
Injection Volume :	0.3	Drying Gas Flow	:5 L/min
Sample Name :	KRAS G12V Peptide (VVGA VGVGK)	T.Flow	:0.2 ml/min
Cat No. :	RP30228	B.conc	:50%H2O/50%MeOH
Lot No. :	E1240424B		
Theoretical MW:	784.95		
Observed MW :	784.6		