

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35
 Cat No.:RP20267
 Lot No. : E1241423J
 Time Processed :8:07:48 AM
 Month-Day-Year Processed :11/02/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.345	18444	2104	0.077
2	8.688	22819219	1240979	95.001
3	9.079	859360	186131	3.578
4	9.575	38667	4005	0.161
5	9.929	46192	3790	0.192
6	10.312	138147	16222	0.575
7	10.588	84423	10362	0.351
8	10.851	15513	2886	0.065
Total		24019967	1466479	100.000