

Sample Name :GRGDSPC

Cat No. : RP20283

Lot No. : E1266826B

Time Processed :15:55:30

Month-Day-Year Processed :02/27/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20510453

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.483	174662	7163	1.068
2	11.600	134045	18255	0.820
3	12.083	15847092	2539755	96.919
4	12.283	125442	20682	0.767
5	12.483	65025	9232	0.398
6	12.825	4644	947	0.028
Total		16350909	2596033	100.000