


Mass Spectrum



Sample Information

Month-Day Processed : 02/28/24
Time Processed : 10:58:16
Injection Volume : 0.3
Sample Name : KRAS WT Peptide (VVVGAGGVGK)
Cat No. : RP30235
Lot No. : E1269724B
Theoretical MW : 842.00
Observed MW : 841.6

Interface : ESI

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

: ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH