

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E1277026B

Time Processed :0:44:47

Month-Day-Year Processed :02/15/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ25010053

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.575	132773	14811	1.374
2	8.767	10467	1823	0.108
3	9.125	9423479	679686	97.534
4	9.558	59516	6990	0.616
5	10.292	20502	1704	0.212
6	17.125	15045	1445	0.156
Total		9661782	706458	100.000