


Mass Spectrum



Sample Information

Month-Day Processed :	02/14/26	Equipment	:ZJ22010150
Time Processed :	20:43:41	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	KRAS G12V Peptide (VVVGAVGVGK)	T.Flow	:0.2 ml/min
Cat No. : RP30229		B.conc	:50%H2O/50%MeOH
Lot No. : E1277026B			
Theoretical MW: 884.08			
Observed MW : 883.8			