


Mass Spectrum



Sample Information

Month-Day Processed :	09/24/22	Equipment	:ZJ22010150
Time Processed :	6:33:56 PM	Interface Bias	:+4.5 kV
Injection Volume :	0.1	Drying Gas Flow	:5 L/min
Sample Name :	β-Endorphin, human	T.Flow	:0.2 ml/min
Cat No. :	RP11344	B.conc	:50%H2O/50%MeOH
Lot No. :	E1349722I		
Theoretical MW :	3465.10		
Observed MW :	3464.8		