


Mass Spectrum



Sample Information

Month-Day Processed :	09/08/22	Equipment	: ZJ21010035
Time Processed :	20:19:17	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	RGD	T.Flow	: 0.2 ml/min
Cat No. :	RP20297	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E13767221		
Theoretical MW:	1025.06		
Observed MW :	1025.0		