

Sample Name :[D-SER]33-Mazdutide
Cat No. : RP30899
Lot No. : E1400525E
Time Processed :9:34:36
Month-Day-Year Processed :06/07/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.833	41747	3506	0.174
2	21.717	856698	52874	3.576
3	22.058	22828079	1351646	95.298
4	22.708	167525	8485	0.699
5	23.175	60314	2899	0.252
Total		23954362	1419409	100.000