

Sample Name : [D-SER]33-Mazdutide
Cat No. : RP30899
Lot No. : E1401425E
Time Processed : 8:23:32
Month-Day-Year Processed : 06/07/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow: 1 ml/min
Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.833	33579	2874	0.145
2	21.725	788789	47122	3.411
3	22.067	22128229	1338648	95.687
4	22.708	125399	6464	0.542
5	23.200	49674	2440	0.215
Total		23125670	1397548	100.000