

Sample Name :[D-GLU]16-Mazdutide
 Cat No. : RP30900
 Lot No. :E1403125E
 Time Processed :17:55:08
 Month-Day-Year Processed :06/10/2025
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ21010375

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.108	40594	2953	0.406
2	19.525	86834	11325	0.869
3	19.783	9688056	903976	96.932
4	20.217	179225	10679	1.793
Total		9994709	928932	100.000