

Sample Name :[D-GLU]16-Mazdutide
Cat No. : RP30900
Lot No. : E1404725E
Time Processed :18:34:30
Month-Day-Year Processed :06/06/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010506

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	19.108	82059	3199	0.870
2	19.425	131231	12236	1.391
3	19.817	9070311	955896	96.119
4	20.608	152931	12059	1.621
Total		9436532	983390	100.000