

Sample Name :[D-ALA]19-Mazdutide
Cat No. : RP30908
Lot No. : E1420825E
Time Processed :15:57:25
Month-Day-Year Processed :06/10/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
10.00	Pumps	B.Conc	35
35.00	Pumps	B.Conc	95
37.00	Pumps	B.Conc	95
37.01	Pumps	B.Conc	5
45.00	Pumps	B.Conc	5
45.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ23010112

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.701	20233	2751	0.288
2	11.832	18287	1888	0.260
3	16.940	6981250	645947	99.275
4	33.733	12429	772	0.177
Total		7032199	651358	100.000