


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	06/06/25	Nebulizing Gas Flow :	1.5L/min	Interface Bias	:ZJ22010150
Time Processed :	6:25:22 PM	CDL Temp :	250	Drying Gas Flow	: +4.5 kV
Injection Volume :	0.6	Block Temp :	200	T.Flow	:5 L/min
Sample Name :	[D-ALA]19-Mazdutide			B.conc	:0.2 ml/min
Cat No. :	RP30908				:50%H2O/50%MeOH
Lot No. :	E1421125E				
Theoretical MW:	4563.14				
Observed MW :	4563.2				