


Mass Spectrum



Sample Information

|                       |                    |                     |           |                 |                 |
|-----------------------|--------------------|---------------------|-----------|-----------------|-----------------|
| Month-Day Processed : | 03/17/23           | Interface           | :ESI      | Equipment       | :ZJ22010150     |
| Time Processed :      | 22:09:16           | Nebulizing Gas Flow | :1.5L/min | Interface Bias  | :+4.5 kV        |
| Injection Volume :    | 0.2                | CDL Temp            | :250      | Drying Gas Flow | :5 L/min        |
| Sample Name :         | Melittin,honey bee | Block Temp          | :200      | T.Flow          | :0.2 ml/min     |
| Cat No. :             | RP20415            |                     |           | B.conc          | :50%H2O/50%MeOH |
| Lot No. :             | E1425423C          |                     |           |                 |                 |

Theoretical MW:2847.50

Observed MW :2847.6