


Mass Spectrum



| | | | | | |
|-----------------------|----------------------|-----------------------|----------|-----------------|-------------------|
| Sample Information | | Interface | | Equipment | |
| Month-Day Processed : | 05/19/25 | Nebulizing Gas Flow : | 1.5L/min | Interface Bias | :ZJ22010150 |
| Time Processed : | 22:25:40 | CDL Temp : | 250 | Drying Gas Flow | : +4.5 kV |
| Injection Volume : | 0.3 | Block Temp : | 200 | T.Flow | :5 L/min |
| Sample Name : | Amylin, human, amide | | | B.conc | :0.2 ml/min |
| Cat No. : | RP11278 | | | | :50% H2O/50% MeOH |
| Lot No. : | E1468125E | | | | |
| Theoretical MW : | 3903.28 | | | | |
| Observed MW : | 3903.3 | | | | |