


Mass Spectrum



Sample Information

|                       |           |                 |                  |
|-----------------------|-----------|-----------------|------------------|
| Month-Day Processed : | 12/26/24  | Equipment       | : ZJ21010035     |
| Time Processed :      | 8:40:25   | Interface Bias  | : +4.5 kV        |
| Injection Volume :    | 0.5       | Drying Gas Flow | : 5 L/min        |
| Sample Name :         | RGD       | T.Flow          | : 0.2 ml/min     |
| Cat No. :             | RP20297   | B.conc          | : 50%H2O/50%MeOH |
| Lot No. :             | E1480524L |                 |                  |
| Theoretical MW:       | 1025.06   |                 |                  |
| Observed MW :         | 1025.0    |                 |                  |