

Sample Name : Glucagon-Like Peptide (GLP) II, human
 Cat No. : RP10769
 Lot No. : E1488622C
 Time Processed : 16:14:39
 Month-Day-Year Processed : 03/20/2022


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	18.434	28949	5297	1.015
2	18.584	9311	1965	0.326
3	18.749	20603	1937	0.722
4	19.097	2719261	314135	95.336
5	19.934	6412	884	0.225
6	22.090	67764	8763	2.376
Total		2852299	332982	100.000