

Sample Name : C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E1539022B
Time Processed : 8:01:26 PM
Month-Day-Year Processed : 02/16/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow: 1 ml/min
Wavelength: 220 nm


Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: GK11010017



1 Detector A Channel 1 / 220nm

Peak Table

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	5.838	9860	25	0.185
2	6.790	24431	1778	0.457
3	7.294	11460	967	0.215
4	7.850	10364	278	0.194
5	8.992	18157	908	0.340
6	12.054	8040	10	0.151
7	12.850	16060	1435	0.301
8	13.671	12015	2951	0.225
9	13.831	5080540	576394	95.104
10	14.129	13031	4039	0.244
11	14.238	4957	798	0.093
12	14.637	103367	16415	1.935
13	18.037	29814	4091	0.558
Total		5342095	610089	100.000