

Sample Name :OVA Peptide (323-339)

Cat No. : RP10610

Lot No. : E1558526C

Time Processed :3:20:10

Month-Day-Year Processed :03/17/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010506

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.625	75614	12100	0.555
2	6.817	26213	4098	0.192
3	7.075	12952478	1246125	95.102
4	7.325	286483	62144	2.103
5	7.467	114013	15988	0.837
6	7.717	102914	11850	0.756
7	7.892	61782	9041	0.454
Total		13619496	1361345	100.000