

Sample Name :Melittin,honey bee
Sample ID :C0769XDCG0-1
Time Processed :18:41:30
Month-Day-Year Processed :08/19/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.00	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010324

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	18.859	6958	738	0.035
2	19.046	18218	2218	0.092
3	19.385	238651	29617	1.210
4	19.752	22212	2112	0.113
5	20.114	225906	24319	1.146
6	20.949	19182474	1873511	97.286
7	21.396	21270	2670	0.108
8	21.858	1961	272	0.010
Total		19717649	1935457	100.000