

Sample Name :Melittin,honey bee
 Cat No. : RP20415
 Lot No. : E1800125C
 Time Processed :21:20:41
 Month-Day-Year Processed :03/19/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.00	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.850	6885	1145	0.090
2	15.317	18117	1492	0.236
3	15.517	31892	1836	0.416
4	15.867	24266	2073	0.317
5	16.167	100357	4915	1.310
6	16.767	7470438	567772	97.479
7	17.267	11704	1887	0.153
Total		7663660	581120	100.000