


Mass Spectrum



Sample Information

Month-Day Processed :	03/20/25	Equipment	:ZJ22010150
Time Processed :	2:47:51 PM	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	Melittin,honey bee	T.Flow	:0.2 ml/min
Cat No. :	RP20415	B.conc	:50%H2O/50%MeOH
Lot No. :	E1800125C		
Theoretical MW:	2847.50		
Observed MW :	2847.6		