


Mass Spectrum



Sample Information

Month-Day Processed : 05/22/24
 Time Processed : 17:39:44
 Injection Volume : 0.3
 Sample Name : DYKDDDDK Peptide
 Cat No. : RP10586
 Lot No. : E1811324E
 Theoretical MW: 1012.97
 Observed MW : 1012.8

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : Z121010035
 Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H2O/50% MeOH