

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)  
Cat No. : RP30229  
Lot No. : E1818825H  
Time Processed :0:52:30  
Month-Day-Year Processed :08/29/2025


Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
Total Flow:1 ml/min  
Wavelength:220 nm

<<LC Time Program>>

| Time  | Module     | Command | Value |
|-------|------------|---------|-------|
| 0.01  | Pumps      | B.Conc  | 5     |
| 25.00 | Pumps      | B.Conc  | 65    |
| 25.01 | Pumps      | B.Conc  | 95    |
| 27.00 | Pumps      | B.Conc  | 95    |
| 27.01 | Pumps      | B.Conc  | 5     |
| 35.00 | Pumps      | B.Conc  | 5     |
| 35.01 | Controller | Stop    |       |

<<Column Performance>>  
<Detector A>  
Column :Inertsil ODS-3 4.6 x 250 mm  
Equipment: ZJ23010113

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area    | Height | Area%   |
|-------|-----------|---------|--------|---------|
| 1     | 9.362     | 28011   | 1828   | 1.155   |
| 2     | 9.615     | 2386177 | 469387 | 98.381  |
| 3     | 12.187    | 11269   | 353    | 0.465   |
| Total |           | 2425457 | 471567 | 100.000 |