

Sample Name:C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E1833024D
Time Processed :3:33:34
Month-Day-Year Processed :04/16/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.225	22615	3138	0.624
2	13.225	16100	2275	0.445
3	13.742	18772	2951	0.518
4	14.650	85780	11754	2.369
5	14.858	34310	3685	0.947
6	15.200	3443885	370569	95.097
Total		3621462	394371	100.000