


Mass Spectrum



Sample Information

Month-Day Processed :	04/15/24	Equipment :	ZJ21010035
Time Processed :	18:27:29	Interface :	ESI
Injection Volume :	0.3	Nebulizing Gas Flow :	1.5L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	CDL Temp :	250
Cat No. :	RP11110	Block Temp :	200
Lot No. :	E1833024D		
Theoretical MW:	2197.60		
Observed MW :	2197.2		

Drying Gas Flow :	5 L/min
T.Flow :	0.2 ml/min
B.conc :	50%H ₂ O/50%MeOH