

Sample Name :TAT peptide
Cat No. :RP20256
Lot No. : E1833125H
Time Processed :17:06:15
Month-Day-Year Processed :08/31/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
31.00	Pumps	B.Conc	95
31.01	Pumps	B.Conc	0
40.00	Pumps	B.Conc	0
40.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010014

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.850	53157	10176	1.971
2	13.000	2571939	524610	95.381
3	13.258	50884	8310	1.887
4	13.992	8585	1837	0.318
5	14.492	5341	1433	0.198
6	14.892	6584	1663	0.244
Total		2696489	548029	100.000