

Sample Name :Influenza A NP(366-374) Strain A/PR/8/35

Cat No. :RP20267

Lot No. :E1835525H

Time Processed :0:40:51

Month-Day-Year Processed :09/02/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010142

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.067	45714	4691	0.535
2	7.858	66257	3836	0.775
3	9.117	36184	4374	0.423
4	9.525	8128580	1181882	95.079
5	9.767	172848	36367	2.022
6	9.983	99743	16877	1.167
Total		8549327	1248026	100.000