


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	06/03/22	Interface	: ZJ21010035
Time Processed :	7:31:46	Nebulizing Gas Flow	: +4.5 kV
Injection Volume :	0.2	CDL Temp	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	Drying Gas Flow	:0.2 ml/min
Cat No. :	RP1111	T.Flow	:50%H2O/50%MeOH
Lot No. :	E1864322E	B.conc	
Theoretical MW :	2197.60		
Observed MW :	2197.2		