

Sample Name :OVA Peptide (323-339)

Cat No. : RP10610

Lot No. : E1882625L

Time Processed :6:21:23

Month-Day-Year Processed :12/23/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ20010139

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.650	4246579	760163	95.152
2	8.842	168795	21513	3.782
3	9.017	35380	6537	0.793
4	11.883	12211	1566	0.274
Total		4462965	789779	100.000